

# AVEA<sup>®</sup> ventilator

## Disposable Expiratory Filter / Watertrap

The HEPA filtration element filters the expired gas to 0.3 microns. With a 99.999+% bacterial and viral efficiency rating, and the heated environment, the expiratory flow sensor minimizes the risk of contamination, while providing an exceptional environment for your staff and patients.

### Resistance comparison

Filter type	Ventilator	Published resistance
Disposable Expiratory Filter / Watertrap (Part number 11790)	AVEA ventilator	< 1.0 cm H <sub>2</sub> O / 60 LPM Flow
Pall with reusable expiratory filter	AVEA ventilator	2.8 cm H <sub>2</sub> O / 100 LPM Flow 1.75 cm H <sub>2</sub> O / 60 LPM Flow

### Benefits of using the AVEA<sup>®</sup> Disposable Expiratory Filter / Watertrap:

- Single patient use
- No component assembly required
- Easy condensation management
- 60% lower resistance than the reusable configuration



Description	Part number	Quantity
AVEA ventilator single patient use expiratory filter / water trap	11790	Case of 12

 **WARNING:** U.S. Federal Law restricts this device to sale by or on the order of a physician.

CareFusion  
22745 Savi Ranch Parkway  
Yorba Linda, CA 92887  
800.231.2466 toll-free  
714.283.2228 tel  
714.283.8493 fax

CareFusion  
Yorba Linda, CA

[carefusion.com](http://carefusion.com)

© 2012 CareFusion Corporation or one of its subsidiaries. All rights reserved. AVEA is a trademark or registered trademark of CareFusion Corporation or one of its subsidiaries. RC3502 (1112/1000) L3479 Rev. A

