

Respiratory ventilation products Critical care

International catalogue



Table of contents

AVEA® ventilator products	1–7
VELA® ventilator products	8–12
Bird® Blender products	13–16
3100A and 3100B High-Frequency Oscillatory Ventilator (HFOV) products	17–22
Infant Flow® nCPAP products.....	23–28
Fisher & Paykel accessories	29–31

AVEA® ventilator products



AVEA accessories



16465
AVEA neonatal hotwire flow sensor
Reusable, each



50000-40035
AVEA water trap for expiratory system
Reusable, each



33985
AVEA water trap bottle for expiratory system
Reusable, each



51000-40640
AVEA removable exhalation cartridge
Reusable, each



33987
AVEA bacterial filter
Reusable, each



22095
AVEA water trap adapter replacement
Reusable, each



51000-02736
Legacy "Bear" circuit support hanger
Each



69300A Circuit support arm with side rail mount bracket (PN 09526)
69300 Circuit support arm without side rail mount bracket (PN 09526)



09526
Side rail mount bracket



34403
Circuit support arm (PN 69300) circuit holder
Pack of 3



15619
Nurse call cable open



15620
Nurse call cable closed



33915
AVEA filter, guard assembly cooling fan
Each



71670
AVEA fan filter
Pack of 5



68289
Sensor, oxygen with connector



11938
Aeroneb nebulizer stand mounting kit
Rail mount bracket (PN 09526) required



Part No.	Description
51000-40897	Fitting, assembly, air inlet, DISS
51000-40888	Fitting, assembly, air inlet, NIST
16144	Fitting, assembly, air inlet, air liquide
16173	Fitting, assembly, air inlet, male DISS
16131	Fitting, assembly, heliox and air, DISS
16132	Fitting, assembly, heliox and air, NIST
16149	Fitting, assembly, heliox and air, air liquide
16174	Fitting, assembly, heliox and air, male DISS

AVEA consumables



51000-40094

AVEA adult flow sensor

Single patient use, each

50000-40031

AVEA adult flow sensor

Single patient use, pack of 10



51000-40098

AVEA neonatal flow sensor

Single patient use, each

50000-40038

AVEA neonatal flow sensor

Single patient use, pack of 10



11395

AVEA exhalation filter


Single patient use, box of 12



11790

AVEA expiratory filter water trap

Single patient use, box of 12




16044

Patient breathing circuit kit, 48"

16045

Patient breathing circuit kit, 72"

AVEA esophageal pressure monitoring consumables




7003100
SmartCath esophageal pressure monitoring tube set
8 Fr, adult
Single patient use, box of 10



7003401
SmartCath esophageal pressure monitoring tube set
6 Fr, pediatric
Single patient use, box of 10



7003300
SmartCath adult nasogastric tube
16 Fr, with esophageal balloon
Single patient use, box of 10



7003402
SmartCath pediatric nasogastric tube
7 Fr, with esophageal balloon
Single patient use, box of 10



79018
Tracheal monitoring tube
Single patient use, each

10635
Tracheal monitoring tube
Single patient use, box of 10



50000-09960
SmartCath extension esophageal pressure monitoring tube
Single patient use, each

50000-09920
SmartCath extension esophageal pressure monitoring tube
Single patient use, box of 10



50000-40034
Tracheal catheter adapter
Single patient use, box of 10

AVEA volumetric capnography accessories



11732

Capnography hardware upgrade kit

Patient cable (PN 16578) and disposable sensors (PN 16605 and 16606) included. Software 4.1 or newer only



16578

Capnostat

CO₂ sensor assembly
Reusable, each



16605

CO₂ adapter, adult/pediatric

Single patient use, box of 10



16606

CO₂ adapter, neonatal

Single patient use, box of 10



16607

CO₂ adapter, adult/pediatric

Reusable, each



16608

CO₂ adapter, neonatal

Reusable, each

CLiO₂ accessories are not sold in the United States



16881
CLiO₂ hardware upgrade kit for CLiO₂-ready AVEA
 SpO₂ module and patient cable (PN 71826) included




To be used with SpO₂ sensor:

71852
LNC-4, LNCS patient cable
 1.2 m (4 ft)
 Each



71853
Infant SpO₂ adhesive sensors
 46 cm (18 in), latex free, weight 3 kg-20 kg
 Single patient use, box of 20



71854
Neonatal SpO₂ adhesive sensors
 46 cm (18 in), latex free,
 weight <3 kg or > 40 kg
 Single patient use, box of 20



71855
Neonatal SpO₂ adhesive sensors
 Latex free, weight < 1 kg
 Single patient use, box of 20




71856
Adult SpO₂ adhesive sensors
 46 cm (18 in), latex free, weight > 30 kg
 Single patient use, box of 20



71857
Pediatric SpO₂ adhesive sensors
 46 cm (18 in), latex free,
 weight 10 kg-50 kg
 Single patient use, box of 20

VELA® ventilator products






16496
VELA Diamond flow sensor
Reusable, each



16240D
VELA Diamond exhalation valve diaphragm
Reusable, pack of 10



10384D
VELA II exhalation valve diaphragm
Reusable, pack of 10



15972
VELA II flow sensor
Reusable, each



20005
Exhalation valve body
(VELA II and Diamond)
Reusable, each



15619
Nurse call cable open



12194
Oximeter mount
For use with transport stand
(PN 11540 or 15987 only)



11539
Exhalation outlet cover kit
For use ONLY with VELA Diamond series,
stands (PN 11536 or PN 11540), and
flow sensor (PN 16496)



15620
Nurse call cable closed

VELA accessories



11938
Aeroneb nebulizer stand mounting kit



69300
Circuit support arm



34403
Circuit support arm (PN 69300) circuit holder
Pack of 3



11536
Basic stand
Support arm not included



21782
Cylinder rings
For use with transport stand (PN 11540)



21773
Basic breakdown stand
Pole, base and casters



21778
Handle
For use with PN 21773



21774
Adapter mounting plate
For use with PN 21773



21775
Adapter mounting post
For use with PN 21773



21780
IV pole



40624
Screw
For use with transport stand (PN 11540)



45038
Split washer
For use with transport stand screw (PN 40624)



16578

Capnostat

CO₂ sensor assembly
Reusable, each



16605

CO₂ adapter, adult/pediatric

Single patient use, box of 10



16607

CO₂ adapter, adult/pediatric

Reusable, each



17509

Partial capnography upgrade kit

Included: power box (PN 16546), reusable adult/pediatric CO₂ adapter (PN 16607), box of 10 single patient use adult/pediatric CO₂ adapter (PN 16605)

Does not include capnography activation software

16920

Full capnography upgrade kit

Included: power box (PN 16546), reusable adult/pediatric CO₂ adapter (PN 16607), box of 10 single patient use adult/pediatric CO₂ adapter (PN 16605)


Includes capnography software with activation



16546

Capnography power box/monitor interface

VELA ventilator circuits



11571

VELA Diamond pediatric breathing circuit kit

Included: circuit tubing, 18" smooth bore; main flow bacteria filter; 90° elbow adapter; 22 mm ID cuff adapter; circuit tubing, 30" smooth bore; watertrap, autoclavable; tapered plug, 7.5 mm male; wye connector; exhalation valve body; Diamond exhalation valve diaphragm

11570

VELA Diamond adult breathing circuit kit

Included: circuit tubing, 18" smooth bore; main flow bacteria filter; 90° elbow adapter; 22 mm ID cuff adapter; circuit tubing, 30" smooth bore; watertrap, autoclavable; tapered plug, 7.5 mm male; wye connector; exhalation valve body; Diamond exhalation valve diaphragm

10684

VELA II adult breathing circuit kit

Included: circuit tubing, 18" smooth bore; main flow bacteria filter; 90° elbow adapter; 22 mm ID cuff adapter; circuit tubing, 30" smooth bore; watertrap, autoclavable; tapered plug, 7.5 mm male; wye connector; exhalation valve body; VELA II exhalation valve diaphragm

Part No.	Description
09531	Circuit tubing, 30" smooth bore (use 09531D for a pack of 10)
09532	Circuit tubing, 18" smooth bore (use 09532D for a pack of 10)
09534	Main flow bacteria filter
04709	90° elbow adapter
09413	Watertrap, autoclavable
4124	Tapered plug, 7.5 mm male
20225	Wye connector
00423	22 mm ID cuff adapter (use 00423D for a pack of 10)

Bird® Blender products



Bird Blenders



03800A
High flow MicroBlender without hoses, DISS

03800B
High flow MicroBlender without hoses, French, Air Liquide inlets

03800J
High flow MicroBlender with hoses, Japan (30 psi differential alarm), DISS

03800G
High flow MicroBlender without hoses, German

10063
High flow MicroBlender without hoses, NIST



10030A
MRI high flow MicroBlender without hoses, DISS



03803T
OmniBlender without hoses, triple outlet, DISS



UltraBlender low flow

Part No.	Description
10040A	0.15–1.0/2–16 LPM, DISS
10041A	0.5–5.0/2–16 LPM, DISS
10042A	1.0–10.0/2–16 LPM, DISS
10046A	2–16/2–16 LPM, DISS
10047A	0.15–1.0/1.0–10.0 LPM, DISS
10048A	0.5–5.0/1.0–10.0 LPM, DISS
10049A	1.0–10.0/1.0–10.0 LPM, DISS
10013A	0.15–1.0/2–16 LPM, NIST
10014A	0.5–5.0/2–16 LPM, NIST
10016A	2–16/2–16 LPM, NIST
10017A	0.15–1.0/1.0–10 LPM, NIST
10018A	0.5–5.0/1.0–10 LPM, NIST



03920A
Low flow MicroBlender without hoses, DISS

10011A
Low flow MicroBlender without hoses, NIST



15625A
Sentry Blender/Analyzer, high flow without hoses, DISS

16694A
Sentry Blender/Analyzer, low flow without hoses, DISS

15625B
Sentry Blender/Analyzer, high flow, French



2688-332
90° elbow adapter



66045
Oxygen sensor, galvanic with flow diverter
For 15 mm port



09220
Flowmeter, 0–15 LPM oxygen
DISS, male



07426
Air inlet water trap/filter assembly



03869H
Auxiliary assembly outlet
9/16–18, 03800



04322
Pole mount blender bracket



05141
Female dovetail bracket



05213
Female dovetail wall bracket



09437
Rail mount blender bracket

Bird Blender accessories



09440
Rail mount bracket



01468
Connector
9/16" (3-way), 2 Male and 1 Female



03837
Connector, tailpiece
Oxygen, 7/16-27



20254
Connector, inlet, air, MicroBlender
NIST



20255
Connector, inlet, oxygen, MicroBlender
NIST



32042
Adapter, oxygen
DISS x 1/4" tube

3100A and 3100B High-Frequency Oscillatory Ventilator (HFOV) products



3100A and 3100B HFOV accessories



467377
Fisher & Paykel humidifier MR850
115 V, with capital orders only

69754
Fisher & Paykel humidifier MR850
230 V, with capital orders only



768965
Humidifier mounting bracket
For use with Teleflex



768968
Humidifier mounting bracket
For use with Fisher & Paykel MR730 and MR850



770134
Fisher & Paykel autofeed pole



770135
Adapter bracket
For use with Fisher & Paykel humidifier base and autofeed pole



770136
Fisher & Paykel
Autofeed pole and bracket kit

Blender and accessories



10273
Blender kit
One blender, air tee/hose assembly and gas output hose



10112
Air tee/hose assembly

Blender and accessories items sold but not shown

Part No.	Description
463605	Hose, oxygen (green)
463606	Hose, air (yellow)
9220	Flow meter 0–15 LPM DISS



766586
Driver cover plate



767138
Support arm for use with adjustable cradle



770566
Circuit cradle




462364
Thumb screw for support arm and cradle

Accessories and service items sold but not shown

Part No.	Description	Qty.
765742	Circuit hold down strap	Each
767163	Air inlet filter element	Pack of 10
766798	Column lint filter element	Each
465298	Krytox inert grease	2 oz tube
765759	3100A HFOV patient circuit calibration label	Each
772754	3100B HFOV patient circuit calibration label	Each

3100A and 3100B HFOV consumables



11437
Connecting tubing set, unfiltered circuits
Red, blue and green tubing
Each



11438
Connecting tubing set, filtered circuits
Red, blue and green tubing
Each



766897
Bellows and water trap set
Box of 4



766896
Cap diaphragm set
Box of 4



12151
Inspiratory filters
Box of 12



12152
Expiratory filters
Box of 25



480029
Circuit port adapter for nitric oxide sampling



772793
Exhaled gas deflector



770150
Scavenger
Export only

Unfiltered circuits



Part No.	Description	Qty.
771374	Fisher & Paykel MR730 humidifier flexible circuit with short heated wire (<i>for heated environment</i>)	Box of 4
771375	Fisher & Paykel MR730 humidifier flexible circuit with long heated wire (<i>for non-heated environment</i>)	Box of 4
773996	Fisher & Paykel MR850 humidifier flexible circuit with short heated wire (<i>for heated environment</i>)	Box of 4
773997	Fisher & Paykel MR850 humidifier flexible circuit with long heated wire (<i>for non-heated environment</i>)	Box of 4
771463	Teleflex 380-90/ConchaTherm III humidifier 16 V, flexible circuit with short heated wire (<i>for heated environment</i>)	Box of 4
771465	Teleflex 380-88/ConchaTherm IV humidifier 21 V, flexible circuit with short heated wire (<i>for heated environment</i>)	Box of 4
771466	Teleflex 380-88/ConchaTherm IV humidifier 21 V, flexible circuit with long heated wire (<i>for non-heated environment</i>)	Box of 4

Unfiltered circuits sold but not shown

Part No.	Description	Qty.
766895	Rigid patient circuit body without heated wire	Box of 4

Filtered circuits



Part No.	Description	Qty.
11518-730S	Fisher & Paykel MR730 humidifier flexible filtered circuit with short heated wire (<i>for heated environment</i>)	Box of 4
11518-730L	Fisher & Paykel MR730 humidifier flexible filtered circuit with long heated wire (<i>for non-heated environment</i>)	Box of 4
11518-850S	Fisher & Paykel MR850 humidifier flexible filtered circuit with short heated wire (<i>for heated environment</i>)	Box of 4
11518-850L	Fisher & Paykel MR850 humidifier flexible filtered circuit with long heated wire (<i>for non-heated environment</i>)	Box of 4
11518-16VS	Teleflex 380-90/ConchaTherm III humidifier 16 V, flexible filtered circuit with short heated wire (<i>for heated environment</i>)	Box of 4
11518-21VS	Teleflex 380-88/ConchaTherm IV humidifier 21 V, flexible filtered circuit with short heated wire (<i>for heated environment</i>)	Box of 4
11518-21VL	Teleflex 380-88/ConchaTherm IV humidifier 21 V, flexible filtered circuit with long heated wire (<i>for non-heated environment</i>)	Box of 4

Notes:

Short heated wires are for heated environments, such as radiant warmers or isolettes
 Long heated wires are for unheated environments, such as cribs, bassinet or beds

Unfiltered circuits

Pre-packaged kit for all required supplies—circuit, bellows, water trap, cap diaphragms and connecting tubing. Box of four complete kits. Humidifier chamber not included.



Part No.	Description	Qty.
11744-730K	Fisher & Paykel MR730 humidifier flexible circuit with heated wire kit	Box of 4
11744-850K	Fisher & Paykel MR850 humidifier flexible circuit with heated wire kit	Box of 4
11744-HRCK	Teleflex 380-88, ConchaTherm IV (21 V), flexible circuit with heated wire, Teleflex ConchaTherm III/IV kit	Box of 4

Filtered circuits

One inspiratory and expiratory filter and a cap diaphragm set. Expiratory filter needs to be changed every 24 hours. Inspiratory filter does not need routine changing.



Part No.	Description	Qty.
16390-101	Fisher & Paykel MR730 humidifier flexible filtered circuit with heated wire	Box of 4
16390-102	Fisher & Paykel MR850 humidifier flexible filtered circuit with heated wire	Box of 4
16390-103	Teleflex 380-88, ConchaTherm IV (21 V) flexible filtered circuit with heated wire	Box of 4

Infant Flow[®] nCPAP products



Infant Flow SiPAP accessories



777237-101

Infant Flow SiPAP stand with IV pole and basket

777237-102

Infant Flow SiPAP stand with IV pole only



777238

IV pole, only

For use with SiPAP stand



777243

Fisher & Paykel MR730 and MR850 humidifier bracket



467377

Fisher & Paykel MR850 humidifier

115 V, with capital orders only

69754

Fisher & Paykel MR850 humidifier

230 V, with capital orders only



466423

Fisher & Paykel MR730 and MR850 humidifier pole clamp



770135

Fisher & Paykel Pole mounting bracket

For attaching auto feed pole



770136

Fisher & Paykel auto feed pole

For use with MR850 humidifier



770134

Fisher & Paykel auto feed pole

For use with MR730 humidifier



466952
Fisher & Paykel MR850 heated wire adapter ("pigtail")

480044
Fisher & Paykel MR559
Heated wire adapter MR730

480038
Fisher & Paykel humidifier adapter kit



467375
Fisher & Paykel MR850 temperature probe
Dual right angle, 72"

467376
Fisher & Paykel MR850 temperature probe
Dual right angle 60"

772211
Fisher & Paykel MR730 temperature probe
Dual right angle, 72"



77248
Teleflex humidifier bracket
Neptune or RCI ConchaTherm humidifier

466448
Teleflex humidifier pole mount clamp
ConchaTherm humidifier



467349
Infant Flow respiratory abdominal sensors
Single patient use, box of 25



677-002
Infant Flow SiPAP transducer assembly



D1420/100
Infant Flow baffle/filter
Single patient use, box of 20



777215
Infant Flow extension tube
Single patient use, pack of 10



775728
Infant Flow bonnet straps
Single patient use, pack of 20


Infant Flow LP nasal CPAP consumables



7772010

Infant Flow LP starter kit

Included: 5 Infant Flow LP generator kits with prongs, 2 XS, S, SM, M, L and XL headgear sizes; 2 XS, S, SM, M, L and XL nasal masks sizes; 2 XL nasal prongs; and 2 XS nasal prongs
Single patient use



7772011

Infant Flow LP generator kit/Fisher & Paykel MR730 or MR850 breathing circuit

Included: Infant Flow LP generator with prongs and neonatal, single-limb, heated breathing circuit

Single patient use, box of 10

7772021

Infant Flow LP generator kit/Fisher & Paykel MR730 or MR850 breathing circuit

Included: Infant Flow LP generator kit with nasal prongs and neonatal, single-limb, heated breathing circuit

Single patient use, box of 20



Nasal masks

For use only with the Infant Flow LP generator

Single patient use, pack of 10

Cat. No.	Color	Size
777002XS	Green	Extra small
777002S	Pink	Small
777002M	Blue	Medium
777002L	Purple	Large
777002XL	Clear	Extra large



Nasal prongs

For use only with the Infant Flow LP generator

Single patient use, pack of 10

Cat. No.	Color	Size
777000XS	Green	Extra small
777000S	Pink	Small
777000M	Blue	Medium
777000L	Purple	Large
777000XL	Clear	Extra large

Headgear

For use with the Infant Flow LP generator assembly
Single patient use



Cat. No.	Size	Measurements	Qty.
777040XS	Extra small	17–21 cm	Box of 10
777040S	Small	21–26 cm	Box of 10
777040SM	Small-medium	24–28 cm	Box of 10
777040M	Medium	26–32 cm	Box of 10
777040L	Large	32–37 cm	Box of 10
777040XL	Extra large	37–42 cm	Box of 10

Bonnets

Single patient use



Cat. No.	Size	Color	Measurements	Qty.
777010	000	White	18–20 cm	Pack of 10
777012	00	Gray	20–22 cm	Pack of 10
777014	0	Pink	22–24 cm	Pack of 10
777016	1	Lt. brown	24–26 cm	Pack of 10
777018	2	Yellow	26–28 cm	Pack of 10
777020	3	Lt. blue	28–30 cm	Pack of 10
777022	4	Gold	30–32 cm	Pack of 10
777024	5	Green	32–34 cm	Pack of 10
777026	6	Lt. burgundy	34–36 cm	Pack of 10
777028	7	Orange	36–38 cm	Pack of 10

Infant Flow nasal CPAP consumables



777085-102

Infant Flow generator kit/circuit

Included: generator, sizing guide and small, medium and large nasal prongs

Single patient use, case of 20

12204-102

Infant Flow generator kit/circuit

Included: generator, prongs and Fisher & Paykel neonatal, single-limb heated circuit for use with MR730 or MR850 humidifier

Single patient use, case of 20

773387-103

Infant Flow generator kit/circuit

Included: generator, prongs and *AirLife* neonatal single-limb heated circuit for use with MR850 humidifier

Single patient use, case of 20

11541-102

Infant Flow generator kit/circuit

Included: generator, prongs and neonatal, single-limb, heated circuit for use with Fisher & Paykel MR730 humidifier

Single patient use, case of 20




Nasal masks

For use only with the Infant Flow generator assembly

Single patient use

Cat. No.	Size	Qty.
777086-101	Small	Box of 10
777086-102	Medium	Box of 10
777086-103	Large	Box of 10
777086-104	Extra large	Box of 10



12204-101

Infant Flow circuit only

Fisher & Paykel neonatal, single-limb heated circuit for use with MR730 or MR850 humidifier

Single patient use, case of 20

773387-101

Infant Flow circuit only

AirLife neonatal, single-limb heated circuit for use with MR850 humidifier and Infant Flow generator

Single patient use, case of 20



Nasal prongs

For use only with the Infant Flow generator assembly

Single patient use

Cat. No.	Size	Qty.
11513-101	Small	Box of 10
11513-102	Medium	Box of 10
11513-103	Large	Box of 10



Bonnets

For use only with the Infant Flow generator assembly

Single patient use, pack of 10

Cat. No.	Size	Color	Measurements
777084-101	000	White	18–20 cm
777084-102	00	Gray	20–22 cm
777084-103	0	Pink	22–24 cm
777084-104	1	Lt. brown	24–26 cm
777084-105	2	Yellow	26–28 cm
777084-106	3	Lt. blue	28–30 cm
777084-107	4	Gold	30–32 cm
777084-108	5	Green	32–34 cm
777084-109	6	Lt. burgundy	34–36 cm
777084-110	7	Orange	36–38 cm
777084-111	8	Dk. green	38–40 cm
777084-112	9	Navy	40–42 cm

Fisher & Paykel accessories





Note: Sold as part of a ventilator sale. Not sold separately.

F&P Healthcare Heated Humidifiers

Part Number	Description
69751	MR850ADU, 230V, Danish
69752	MR850AEA, 230V, English
69753	MR850AEK, 230V, English
69754	MR850AEU, 230V, English/Spanish/French
69755	MR850AFU, 230V, French/Dutch/German
69756	MR850AGU, 230V, German/French/Italian
69757	MR850ALU230V, Italian/Spanish/Portuguese
69758	MR850ANU, 230V, Norwegian/Swedish/Finnish
69759	MR850ARU, 230V, English/Russian/Polish
69760	MR850GJU, 100V, Japanese
69761	MR850JHU, 115V, English/Spanish/French
69762	MR850JSU, 115V, Spanish/Portuguese/English



F&P Humidification Chambers

Part Number	Description
69323	MR225 neonatal SPU box of 50
69319	MR290U autofeed SPU box of 40
69318	MR290UX autofeed SPU box of 10
69370	MR370 adult, reuse



Humidifier Airway Temperature Probes

Part Number	Description
69308	900MR561 dual airway temp probes, 2.2 m (88") MR700 series/HC500
69764	900MR869 dual airway temp probes, 1.85 m (74") MR800 series
69769	900MR861 dual airway temp probes, 2.2m (88")

Heater wire adapters



69763
900MR805 electrical adapter, dual limb Fisher & Paykel circuit, MR800 series



466952
900MR806 electrical adapter, single limb Fisher & Paykel circuit, MR800 series



69558
900MR558 electrical adapter, single limb Lemo style for reusable circuit, MR700 series



69765
900MR858 electrical adapter, dual limb Lemo style for reusable circuit, MR700 series



69301
Bracket, Bird parallel channel adapter



466423
Pole mount C-clamp 19 to 40 mm pole




69302
Pole mount bracket 15 to 50 mm



770135
Bracket adapter for mounting pole



69316
Heater wire assembly, 1.1 m (43"), bias flow > 4 lpm



Draw wire (85"), reuse 69310

Part Number	Description
69309	Draw wire (60"), reuse
69310	Draw wire (85"), reuse



69306
Adapter, 10 mm M to 15 mm M / 10 mm F reuse



69131
Elbow adapter with temp probe port, 22 mm F to 22 mm M / 15 mm F reuse

Items sold, but not shown

Part No.	Description
770134	Water bag pole, bent, 1.05 m (42")
69515	Heater wire assembly, 1.5 m (60")



CareFusion
Yorba Linda, CA

carefusion.com/criticalcare



© 2012 CareFusion Corporation or one of its subsidiaries. AVEA, VELA, EnVe, Bird Blenders, Infant Flow, AirLife, CLIO₂, CareFusion and the CareFusion logo are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries. RC4090 (1012/PDF) L3508